



THE UNIVERSITY of EDINBURGH
Royal (Dick) School of
Veterinary Studies



Dentistry Equipment and Instrumentation: Organising Your Dentistry Workspace



Claire Harrison BSc (Hons) RVN VTS (Dentistry) AFHEA
Royal (Dick) School of Veterinary Studies, University of Edinburgh & iM3
Dental Ltd; Email: Claire.Harrison_vtsdent@yahoo.com

Dental Equipment: Why is it so important?

Dentistry caseload in practice- Prevalence of dental conditions. 80% dogs & 70% cats over 3 years old will have periodontal disease.

Right tools for the right job!

Improve efficiency and standard of dentistry and less operator fatigue.

AAHA Dental Guidelines 2019. www.aaha.org/globalassets/02-guidelines/dental/aaha_dental_guidelines.pdf



The Dental Operating Space

- Common limitations- Space & location, equipment position, lighting, table design.
- Consider ergonomic design
- Smoother, more effective procedures.
- Less damage to equipment and instruments.



Ergonomics in the Dental Operating Space

- Definition of Ergonomics- (Human Factors)
- Application of psychological and physiological principle during design of products, processes and systems.
- Goal- enhance safety and comfort and interaction, increase productivity. Reduce human error.





Improving Ergonomics




- Table design
- Equipment location
- Patient positioning
- Lighting & Magnification

The Modified Pen-Grip.



Palm Grip



Personal Protective Equipment (PPE)

- Gloves
- Mask
- Goggles/ Eyeshield
- Hair net
- Theatre wear



Patient Comfort & Safety

- Spring loaded gags push TMJ's beyond ROM
- Metal gags damage tooth enamel
- Consider soft rubber wedges
- Relieve gag pressure every 5 mins
- Cats!!



Dental Machines- How do They Work?

- ▶ Air driven compressor.
- ▶ Oil Vs Oil Free compressor



Dental Machines- Maintenance & Infection Control

- Daily/weekly- Clean after use, change distilled water & release compressor
- Annual Service- qualified engineer
- Infection control- biofilm- What can we do?
- Troubleshooting & Create Dental SOP's



Ultrasonic Scalers



- Magnetostrictive- move in predictable pattern. Elliptical
- Magnetic field energises water flow on scaler tip.
- Higher frequency range – 40,000Hz.
- Piezoelectric- Linear movement. Adjustable frequency.



Identifying- Dental Handpieces




Handpiece care, oiling and sterilisation protocol.

!!No Cross Contamination!!




Oiling protocol for Handpieces





Equipment for the Oral Assessment.



- Separating and organising instruments for different stages of the procedure.
- What do we need available?

FELINE DENTAL CHART

Pfizer Animal Health

DATE	
ANIMAL'S NAME	AGE
OWNER'S NAME	SEX
BREED	NUMBER

Buccal

R

MAXILLA

L

Palatal

109 108 107 106 104 103 102 101 201 202 203 204 206 207 208 209

CI															
GI															
Other															

Buccal

R

MANDIBLE

L

Lingual

409 408 407 404 403 402 401 301 302 303 304 207 208 209

CI															
GI															
Other															

Treatment	Recommendations
-----------	-----------------

GI (gingivitis index)

0 = Normal gingiva

1 = Marginal gingivitis, not bleedable

2 = Bleeds on gentle probing, swollen

3 = Severe inflammation, spontaneous bleeding

CI (calculus index)

0 = No calculus

1 = Up to 25% crown buccally/lingually

2 = From 25% to 75% crown on buccal/lingual

3 = From 75% to 100% crown on buccal/lingual

Feline's FICIR - stain type and class eg FICIR 2/4 Requires X-ray

Type 1: Local infection in crown/roots, class of root/s 1-4

Type 2: Local infection of roots & replacement by bone

Class 1 = Crown intact

Class 2 = Penetration into dentine

Class 3 = Penetration into pulp

Class 4 = Major destruction of root

Class 5 = Crown fractured/severed

X-ray of tooth taken:

AX - Alveolar

CM - Crown

OD - Ocular

GWP - Gingivitis

CODE KEY

IR - Missing tooth

ER - Extracted tooth

PR - Preserved tooth

FX - Fractured tooth

F - Feline's FICIR

GI - Gingivitis

OD - Ocular

CM - Crown

AX - Alveolar

GWP - Gingivitis

Published by Pfizer Animal Health with the kind assistance of Norman W. Johnson, WSAVA World Dental Society member

DENTAL LEADERSHIP



Dental Radiography



Periodontal Therapy Kit-What is it for?



- Diagnostic
- Periodontal cleaning and treatment.
- What sterilised instruments should we have?

Periodontal Probes

- Blunt round ended instrument.
- Indentations/coloured bands for measurements.
- Used for gingivitis score and measurements of periodontal disease.
- Sizes and types- 10, 12 and 16mm



Dental explorer

- Sharp pointed tip.
- For use on tooth surface only.
- Find defects, Caries and tooth fractures.
- Types e.g. Shepherds Hook



Dental Mirror and Handle



Supragingival Scaler

- For use on crown of tooth only
- Keep away from gingiva
- Sharp working edges
- Cleaning developmental grooves and fissures and interdental spaces.



Subgingival Curettes




- Used for subgingival curettage and root planing.
- Used under the gumline i.e perio pockets.
- Curved rounded back
- Gracey curettes- one working edge.

Surgical Extraction Kit



- **Ophthalmic instruments**
- **High quality**
- **In tray autoclaved & ready**





Dental Extraction & Oral Surgery Instruments



Instrument selection – Extractions

- **Instruments to make and close flaps**
- **Burs to facilitate extraction**
- **Luxators to cut PL fibres and loosen roots**
- **Selection of luxator blade sizes**
- **Extraction forceps**
- **All sharpened, sterilised & ready to use**




Periosteal Elevators

**Molt 2/4
periosteal
elevator**



Surgical Scissors



- Very sharp- separating periosteum.
- Types- Iris, Goldman fox, La Grange
- Metzenbaum- blunt dissection
- Dissecting Forceps.

Dental Highspeed Burs



- Friction Grip (FG)
- Tungsten Carbide or Stainless steel.
- Cross cut Taper Fissure-701, 701L
- Round Bur- size 1,2,3 &4



Diamond Burs



- Friction grip (FG)
- Various sizes & shapes,
- Round & rugby ball shaped most useful.
- Grit- Fine - course (recommend medium)
- Used for- alveoloplasty, odontoplasty.

Luxator or Elevator?



Elevator

**Thick blade
Used like a shoehorn
Mostly kept above
bone crest**


Luxator

**Thin sharp blades to
insert into PL space.
Must be sharp &
well maintained**

Feline Extractions



- Slimline luxating elevators.
- Modified versions
- Atraumatic root pick
- Winged elevators



Extraction Forceps



- Recommend small breed sized forceps.
- Fully closing.
- Cats & Root tips?


Instrument Sharpening



- Sharpening oil or water – dependent on stone type.
- Stationary stone-moving instrument technique.
- Ensure correct angle for instrument. Gentle pressure required.
- Sharpening Machine!

When and What to Sharpen

- When- After every use.
Results in smoother procedures & increased instrument life.
- Hand scalers & Curettes.
Periosteal elevators &
Luxators/elevators.



Vet-Tome



Equipment advances in Dentistry & Oral Surgery.



- Piezoelectric devices- Ultrasonic bone cutting instruments.
- High powered, low frequency= specific to cutting mineralised tissue.
- Reduced trauma to soft tissue and risk of microfractures compared to oscillating saws and drills.

Alternative Diagnostic Imaging



- Handheld Radiography Units.
- Cone Beam Computed Tomography (CBCT). Increased diagnostic yield for dentoalveolar pathology.
- Limitations

- Ability to evaluate the dental operating space and necessary requirements to improve ergonomics.
- Understand the importance of dental equipment and instrumentation in effective dental procedures.
- Have the ability to identify commonly used dental equipment and instruments. It's care, organisation and maintenance of said equipment.
- Understand the principles behind correct sharpening techniques.



References & Further Reading

- Angel M, 2014, Sharpening Periodontal Instruments, J Vet Dent, Vol:31(1), pg 58-64.
- Ed Reiter AM, Gracis M, BSAVA Manual of Canine and Feline Dentistry and Oral Surgery 4th Ed, 2018, BSAVA, Gloucester.
- Perrone, JR, Small Animal Dental Procedures for Veterinary Technicians & Nurses, 2013, Wiley-Blackwell, Singapore.

